Professional Instruction

Kinesiologic Tape cervical spine

1. Application

Thanks to their adhesive wick is applied in a wavy pattern, the tapes adhere pleasantly to the skin. The tapes can and should be worn for several days. After approximately 5-7 days, their position should be examined and new tapes should be reapplied if necessary. More extended wear has no harmful consequences, even without examination.

Verifying that the skin is sound and unbroken is recommended prior to each application.

2. Areas of application

Area of the upper spine (cervical spine). In case of radiating pain, especially in combination with sensitivity or motor disturbances, consulting a doctor or physiotherapist is highly recommended.

3. Precautions, safety information

When the complaint is first noticed, the cause should be examined medically in order to exclude the possibility of causes that require medical treatment.

The attending doctor or a physiotherapist should be consulted in case of severe or long-term complaints.

4. Particularities of the tapes

The tapes do not contain harmful additives. The adhesive that is used has been dermatologist tested and is considered especially skin-friendly. This is why there is also normally no risk for persons allergic to conventional plasters. Combinations with other types of therapy are possible and/or recommended (supporting therapy).

Thanks to the high level of wearer comfort of the tapes, all day-to-day work and sports activities can and/or should be conducted as usual with the tapes in place. The same applies to showering. The material is breathable, permeable to moisture and liquids and dries quickly.

5. Application method and removal of the tapes

The skin must be dry and free of grease. Hair should be trimmed (long-hair trimmer) or gently removed entirely by shaving (if possible, the day before applying the tapes). Intensively pressing the tapes in place and rubbing (especially along the edges) after application ensures the adhesive will adhere properly.

To remove the tapes, carefully and gradually pull in the direction of hair growth. The tapes are especially easy to remove in the shower or while taking a bath. Skin cream or body oil help relax the skin after removing the tapes.

The cervical spine tapes consist of two narrow strips, approximately 30 cm in length. However, the length is individually adjustable. After adjusting the length, the edges have to be rounded again.

Both strips are applied to the shoulder-neck region.

Take the strips out of the package and stand upright.

1. Remove the backing paper from the first tape strip. Avoid touching the adhesive surface during this process.

2. Apply one end of a tape strip to the bony edge of the right shoulder (acromion) so it is centred and secure this starting point.

3. Now push the cervical spine region out slightly by pressing your chin to your chest and slightly tilting your head towards the left shoulder or looking at the tape starting point until you feel slight tension or you reach the limit of your range of motion. If you experience discomfort or the level of existing discomfort increases, stop the movement and raise the head again until the level of discomfort decreases or the feeling goes away.

4. Now apply the tape without tension in the direction of the occiput (back of the head) and press it into place. Make sure no folds or ripples form in the process.

5. Proceed the same way with the second tape strip on the other shoulder / side of the neck.

6. Precautions, safety information

- » Always cut the tapes to size prior to application, never while you are applying them.
- » Always be careful when pulling clothing over the edges of the tapes; avoid loosening the tapes.
- » Avoid the formation of folds and ripples, especially where they could cause pressure.
- » If your complaints persist or worsen, consult your doctor and / or physiotherapist. Never apply the tapes if your complaints have not been diagnosed.
- » Remove the tapes immediately in case of intolerance.
- » Store the material at room temperature and avoid direct exposure to sunlight.

7. Side effects

There are no known side effects with correct application.

BIOFREEZE

Biofreeze provides temporary relief from many aches and pains, such as sore muscles and joints, sprains, strains, and bruises. It can be applied to body parts such as the foot, ankle, knee, hand, wrist, elbow, hip, back, neck and shoulders. The Biofreeze Original Gel is particularly well-suited for massage applications on larger muscle areas that may benefit from a hands-on approach.

- Penetrates quickly to relieve aches and pains
- Proprietary formula made with USP-grade ingredients
- Free of propylene glycol, dyes and parabens

For external use only. Do not apply to wounds or damaged skin. Avoid contact with the eyes or mucous membranes.

Note:

Treatment with kinesiology tape is based on methods from training theory and kinesiology studies, the effects of which are not proven by means of conventional medical studies but rather are assumed based on experience reports that are not evidence-based as well as application observations.

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© 2018 Ludwig Artzt GmbH.
Ludwig Artzt GmbH
© 2017 Ludwig Artzt GmbH.
©2020 Ludwig Artzt GmbH. Alle Rechte geschützt.